




DENND2C Antibody

Product Code	CSB-PA724250LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68D51
Immunogen	Recombinant Human DENN domain-containing protein 2C protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DENN domain-containing protein 2C, DENND2C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DENND2C

Image



Western blot

All lanes: DENND2C antibody at 8µg/ml + A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 107, 95, 101 kDa

Observed band size: 107 kDa