




LGALS9C Antibody

Product Code	CSB-PA728553LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6DKI2
Immunogen	Recombinant Human Galectin-9C protein (1-356AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-9C (Gal-9C) (Galectin-9-like protein B), LGALS9C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LGALS9C

Image



Western blot
 All lanes: LGALS9C antibody at 4µg/ml
 Lane 1: Mouse stomach tissue
 Lane 2: Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 40 kDa
 Observed band size: 22 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.