




RTP2 Antibody

Product Code	CSB-PA728807LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5QGT7
Immunogen	Recombinant Human Receptor-transporting protein 2 protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-transporting protein 2 (3CxxC-type zinc finger protein 2), RTP2, Z3CXXC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RTP2

Image



Western blot
 All lanes: RTP2 antibody at 8µg/ml + Mouse liver tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 27 kDa
 Observed band size: 27 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.