



L1TD1 Antibody, Biotin conjugated

Product Code	CSB-PA729208LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T7N2
Immunogen	Recombinant Human LINE-1 type transposase domain-containing protein 1 protein (665-865AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LINE-1 type transposase domain-containing protein 1 (ES cell-associated protein 11), L1TD1, ECAT11
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	L1TD1