




DMRTA1 Antibody

Product Code	CSB-PA730143LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VZB9
Immunogen	Recombinant Human Doublesex- and mab-3-related transcription factor A1 protein (201-504AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Doublesex- and mab-3-related transcription factor A1, DMRTA1, DMO
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DMRTA1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA730143LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.