



Lcn3 Antibody, Biotin conjugated

Product Code	CSB-PA730794HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62471
Immunogen	Recombinant Mouse Vomeronasal secretory protein 1 protein (19-182AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vomeronasal secretory protein 1 (Lipocalin-3) (Vomeronasal secretory protein I) (VNSP I), Lcn3
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Lcn3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.