



Lcn3 Antibody, Biotin conjugated

Product Code	CSB-PA730794HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62471
Immunogen	Recombinant Mouse Vomeronasal secretory protein 1 protein (19-182AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vomeronasal secretory protein 1 (Lipocalin-3) (Vomeronasal secretory protein I) (VNSP I), Lcn3
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Lcn3