



ZKSCAN2 Antibody, Biotin conjugated

Product Code	CSB-PA730959LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63HK3
Immunogen	Recombinant Human Zinc finger protein with KRAB and SCAN domains 2 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 2 (Zinc finger protein 694), ZKSCAN2, ZNF694
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZKSCAN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.