



Cyp3a13 Antibody

Product Code	CSB-PA731056LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q64464
Immunogen	Recombinant Mouse Cytochrome P450 3A13 protein (1-503AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 3A13 (EC 1.14.14.1) (CYPIIIA13), Cyp3a13, Cyp3a-13
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cyp3a13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.