





TEX35 Antibody

Product Code	CSB-PA732917LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T0J7
Immunogen	Recombinant Human Testis-expressed protein 35 protein (1-215AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis-expressed protein 35, TEX35, C1orf49
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TEX35

Image



Western blot
All lanes: TEX35 antibody at 4µg/ml + U251 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 27, 26, 24, 21, 25 kDa
Observed band size: 27 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA732917LA01HU at dilution of 1:100

