





DDI2 Antibody

Product Code	CSB-PA732937LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5TDH0
Immunogen	Recombinant Human Protein DDI1 homolog 2 protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein DDI1 homolog 2 (EC 3.4.23.-), DDI2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DDI2

Image



Western blot
 All lanes: DDI2 antibody at 4µg/ml
 Lane 1: A431 whole cell lysate
 Lane 2: HL60 whole cell lysate
 Lane 3: K562 whole cell lysate
 Lane 4: Caco-2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 45, 24, 47 kDa
 Observed band size: 47 kDa



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA732937LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.