



OLFML2B Antibody, HRP conjugated

Product Code	CSB-PA734650LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68BL8
Immunogen	Recombinant Human Olfactomedin-like protein 2B protein (451-750AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactomedin-like protein 2B (Photomedin-2), OLFML2B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OLFML2B