



NCMAP Antibody, Biotin conjugated

Product Code	CSB-PA735965LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T1S8
Immunogen	Recombinant Human Noncompact myelin-associated protein (54-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Noncompact myelin-associated protein (Myelin protein of 11 kDa) (MP11), NCMAP, C1orf130
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NCMAP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.