



# NCMAP Antibody, Biotin conjugated

Product Code	CSB-PA735965LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T1S8
Immunogen	Recombinant Human Noncompact myelin-associated protein (54-102AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Noncompact myelin-associated protein (Myelin protein of 11 kDa) (MP11), NCMAP, C1orf130
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NCMAP