



FAM102B Antibody, FITC conjugated

Product Code	CSB-PA735976LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T8I3
Immunogen	Recombinant Human Protein FAM102B protein (92-360AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM102B, FAM102B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM102B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.