



MORN5 Antibody, Biotin conjugated

Product Code	CSB-PA736208LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VZ52
Immunogen	Recombinant Human MORN repeat-containing protein 5 protein (1-161AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORN repeat-containing protein 5, MORN5, C9orf113 C9orf18
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MORN5