



ATG1 Antibody, HRP conjugated

Product Code	CSB-PA738902LB01CZI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6FL58
Immunogen	Recombinant Candida glabrata Serine/threonine-protein kinase ATG1 protein (11-312AA)
Raised In	Rabbit
Species Reactivity	Candida glabrata
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase ATG1 (EC 2.7.11.1) (Autophagy-related protein 1), ATG1
Immunogen Species	Candida glabrata
Research Area	Others
Target Names	ATG1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.