



ATG8 Antibody, Biotin conjugated

Product Code	CSB-PA739788HD01CZI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6FXR8
Immunogen	Recombinant Candida glabrata Autophagy-related protein 8 protein (1-116AA)
Raised In	Rabbit
Species Reactivity	Candida glabrata
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Autophagy-related protein 8 (Autophagy-related ubiquitin-like modifier ATG8), ATG8, AUT7
Immunogen Species	Candida glabrata
Research Area	Others
Target Names	ATG8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.