



acpS Antibody, HRP conjugated

Product Code	CSB-PA740005HB01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GF02
Immunogen	Recombinant Staphylococcus aureus Holo-[acyl-carrier-protein] synthase protein (1-119AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Holo-[acyl-carrier-protein] synthase (Holo-ACP synthase) (EC 2.7.8.7) (4'-phosphopantetheinyl transferase AcpS), acpS, dpj
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Research Area	Others
Target Names	acpS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.