



TBC1D3C Antibody, HRP conjugated

Product Code	CSB-PA740248LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IPX1
Immunogen	Recombinant Human TBC1 domain family member 3C protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 3C, TBC1D3C, TBC1D3D
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TBC1D3C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.