



# OPN5 Antibody, Biotin conjugated

Product Code	CSB-PA740895LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6U736
Immunogen	Recombinant Human Opsin-5 protein (310-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opsin-5 (G-protein coupled receptor 136) (G-protein coupled receptor PGR12) (Neuropsin) (Transmembrane protein 13), OPN5, GPR136 PGR12 TMEM13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPN5