



OPN5 Antibody, Biotin conjugated

Product Code	CSB-PA740895LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6U736
Immunogen	Recombinant Human Opsin-5 protein (310-353AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Opsin-5 (G-protein coupled receptor 136) (G-protein coupled receptor PGR12) (Neuroopsin) (Transmembrane protein 13), OPN5, GPR136 PGR12 TMEM13
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OPN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.