



ENPP6 Antibody, HRP conjugated

Product Code	CSB-PA740915LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UWR7
Immunogen	Recombinant Human Ectonucleotide pyrophosphatase/phosphodiesterase family member 6 protein (31-330AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ectonucleotide pyrophosphatase/phosphodiesterase family member 6 (E-NPP 6) (NPP-6) (EC 3.1.4.-) (EC 3.1.4.38) (Choline-specific glycerophosphodiester phosphodiesterase) (Glycerophosphocholine cholinephosphodiesterase) (GPC-Cpde), ENPP6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ENPP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.