



FAM151B Antibody, Biotin conjugated

Product Code	CSB-PA740923LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXP7
Immunogen	Recombinant Human Protein FAM151B protein (1-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM151B, FAM151B
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM151B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.