



ZNF767P Antibody, Biotin conjugated

Product Code	CSB-PA741699LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q75MW2
Immunogen	Recombinant Human Protein ZNF767 protein (1-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ZNF767 (Zinc finger protein 767 pseudogene), ZNF767P, ZNF767
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF767P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.