



KCTD9 Antibody, Biotin conjugated

Product Code	CSB-PA741970LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L273
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD9 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing protein KCTD9, KCTD9
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCTD9