



GALNT17 Antibody, HRP conjugated

Product Code	CSB-PA743592LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IS24
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 17 protein (101-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polypeptide N-acetylgalactosaminyltransferase 17 (EC 2.4.1.41) (Polypeptide GalNAc transferase-like protein 3) (GalNAc-T-like protein 3) (pp-GaNTase-like protein 3) (Protein-UDP acetylgalactosaminyltransferase-like protein 3) (UDP-GalNAc:polypeptide N-acetylgalactosaminyltransferase-like protein 3) (Williams-Beuren syndrome chromosomal region 17 protein), GALNT17, WBSCR17
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GALNT17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.