




SLC30A6 Antibody

Product Code	CSB-PA744010LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NXT4
Immunogen	Recombinant Human Zinc transporter 6 protein (246-461AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 6 (ZnT-6) (Solute carrier family 30 member 6), SLC30A6, ZNT6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC30A6

Image



Western blot

All lanes: SLC30A6 antibody at 6µg/ml

Lane 1: A549 whole cell lysate


Lane 2: U251 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 52, 56, 48, 49 kDa


Observed band size: 52 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA744010LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA744010LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA744010LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.