



SLC30A6 Antibody, Biotin conjugated

Product Code	CSB-PA744010LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NXT4
Immunogen	Recombinant Human Zinc transporter 6 protein (246-461AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter 6 (ZnT-6) (Solute carrier family 30 member 6), SLC30A6, ZNT6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC30A6