




SFXN4 Antibody

Product Code	CSB-PA744046LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4A7
Immunogen	Recombinant Human Sideroflexin-4 protein (2-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC, IF; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sideroflexin-4 (Breast cancer resistance marker 1), SFXN4, BCRM1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SFXN4

Image



Western blot

All lanes: SFXN4 antibody at 0.3μg/ml

Lane 1: K562 whole cell lysate

Lane 2: U251 whole cell lysate

Lane 3: U87 whole cell lysate

Lane 4: A549 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 38, 35, 25 kDa


Observed band size: 38 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA744046LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA744046LA01HU at dilution of 1:100



Immunofluorescent analysis of U251 cells using CSB-PA744046LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.