



SFXN4 Antibody, Biotin conjugated

Product Code	CSB-PA744046LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4A7
Immunogen	Recombinant Human Sideroflexin-4 protein (2-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sideroflexin-4 (Breast cancer resistance marker 1), SFXN4, BCRM1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SFXN4