



SFXN4 Antibody, Biotin conjugated

Product Code	CSB-PA744046LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4A7
Immunogen	Recombinant Human Sideroflexin-4 protein (2-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sideroflexin-4 (Breast cancer resistance marker 1), SFXN4, BCRM1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SFXN4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.