



KCTD1 Antibody, Biotin conjugated

Product Code	CSB-PA744477LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q719H9
Immunogen	Recombinant Human BTB/POZ domain-containing protein KCTD1 protein (1-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BTB/POZ domain-containing protein KCTD1 (Potassium channel tetramerization domain-containing protein 1), KCTD1, C18orf5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	KCTD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.