



PF14_0076 Antibody, Biotin conjugated

Product Code	CSB-PA745322LD01EWP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7KQM4
Immunogen	Recombinant Plasmodium falciparum Plasmeprin-1 protein (125-452AA)
Raised In	Rabbit
Species Reactivity	Plasmodium falciparum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasmeprin-1 (EC 3.4.23.38) (Aspartic hemoglobinase I) (PfAPG) PF14_0076
Immunogen Species	Plasmodium falciparum
Research Area	Others
Target Names	PF14_0076
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.