



THOC7 Antibody, FITC conjugated

Product Code	CSB-PA746911LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6I9Y2
Immunogen	Recombinant Human THO complex subunit 7 homolog protein (4-204AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	THO complex subunit 7 homolog (Functional spliceosome-associated protein 24) (fSAP24) (Ngg1-interacting factor 3-like protein 1-binding protein 1) (NIF3L1-binding protein 1) (hTREX30), THOC7, NIF3L1BP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	THOC7