



THOC7 Antibody, Biotin conjugated

Product Code	CSB-PA746911LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6I9Y2
Immunogen	Recombinant Human THO complex subunit 7 homolog protein (4-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	THO complex subunit 7 homolog (Functional spliceosome-associated protein 24) (fSAP24) (Ngg1-interacting factor 3-like protein 1-binding protein 1) (NIF3L1-binding protein 1) (hTREX30), THOC7, NIF3L1BP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	THOC7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.