




PYHIN1 Antibody

Product Code	CSB-PA746957LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6K0P9
Immunogen	Recombinant Human Pyrin and HIN domain-containing protein 1 protein (1-245AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyrin and HIN domain-containing protein 1 (Interferon-inducible protein X), PYHIN1, IFIX
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PYHIN1

Image



Western blot

All lanes: PYHIN1 antibody at 6µg/ml + Mouse small intestine tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 56, 55, 52, 51, 27, 26 kDa

Observed band size: 56, 22, 84 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.