



RINT1 Antibody, Biotin conjugated

Product Code	CSB-PA747351LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUQ1
Immunogen	Recombinant Human RAD50-interacting protein 1 protein (401-700AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAD50-interacting protein 1 (RAD50 interactor 1) (HsRINT-1) (RINT-1), RINT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RINT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.