




MORN3 Antibody

Product Code	CSB-PA747476LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PF18
Immunogen	Recombinant Human MORN repeat-containing protein 3 protein (1-240AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORN repeat-containing protein 3, MORN3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MORN3

Image



Western blot
 All lanes: MORN3 antibody at 2µg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 28, 19 kDa
 Observed band size: 28 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.