



WDR74 Antibody, Biotin conjugated

Product Code	CSB-PA747557LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6RFH5
Immunogen	Recombinant Human WD repeat-containing protein 74 protein (111-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 74 (NOP seven-associated protein 1), WDR74, NSA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	WDR74
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.