



HTR3D Antibody, Biotin conjugated

Product Code	CSB-PA747864LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70Z44
Immunogen	Recombinant Human 5-hydroxytryptamine receptor 3D protein (107-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	5-hydroxytryptamine receptor 3D (5-HT3-D) (5-HT3D) (Serotonin receptor 3D), HTR3D
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HTR3D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.