




# fabI Antibody

<b>Product Code</b>	CSB-PA750140HA01SKV
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q6GI75
<b>Immunogen</b>	Recombinant Staphylococcus aureus Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI protein (1-256AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Staphylococcus aureus, Mouse
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI (ENR) (EC 1.3.1.39) (NADPH-dependent enoyl-ACP reductase), fabI
<b>Immunogen Species</b>	Staphylococcus aureus
<b>Research Area</b>	Others
<b>Target Names</b>	fabI

## Image



Western blot  
 All lanes: fabI antibody at 5µg/ml + Mouse heart tissue  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 28 kDa  
 Observed band size: 28 kDa

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.