




fabI Antibody

Product Code	CSB-PA750140HA01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GI75
Immunogen	Recombinant Staphylococcus aureus Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI protein (1-256AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI (ENR) (EC 1.3.1.39) (NADPH-dependent enoyl-ACP reductase), fabI
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	fabI

Image



Western blot
All lanes: fabI antibody at 5µg/ml + Mouse heart tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 28 kDa
Observed band size: 28 kDa