



fabI Antibody, HRP conjugated

Product Code	CSB-PA750140HB01SKV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6GI75
Immunogen	Recombinant Staphylococcus aureus Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI protein (1-256AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Enoyl-[acyl-carrier-protein] reductase [NADPH] FabI (ENR) (EC 1.3.1.39) (NADPH-dependent enoyl-ACP reductase), fabI
Immunogen Species	Staphylococcus aureus (strain MRSA252)
Research Area	Others
Target Names	fabI
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.