



SNRNP48 Antibody

Product Code	CSB-PA750321LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6IEG0
Immunogen	Recombinant Human U11/U12 small nuclear ribonucleoprotein 48 kDa protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U11/U12 small nuclear ribonucleoprotein 48 kDa protein (U11/U12 snRNP 48 kDa protein) (U11/U12-48K), SNRNP48, C6orf151
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SNRNP48
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.