





GRHL2 Antibody

Product Code	CSB-PA750350LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ISB3
Immunogen	Recombinant Human Grainyhead-like protein 2 homolog protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Grainyhead-like protein 2 homolog (Brother of mammalian grainyhead) (Transcription factor CP2-like 3), GRHL2, BOM TFPCP2L3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRHL2


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA750350LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA750350LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA750350LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.