



GRHL2 Antibody, Biotin conjugated

Product Code	CSB-PA750350LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ISB3
Immunogen	Recombinant Human Grainyhead-like protein 2 homolog protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Grainyhead-like protein 2 homolog (Brother of mammalian grainyhead) (Transcription factor CP2-like 3), GRHL2, BOM TF2L3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRHL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.