



FAM71F2 Antibody, Biotin conjugated

Product Code	CSB-PA750765LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NXP2
Immunogen	Recombinant Human Protein FAM71F2 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM71F2 (Protein FAM137B), FAM71F2, FAM137B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM71F2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.