



ARHGAP11A Antibody

Product Code	CSB-PA750800LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P4F7
Immunogen	Recombinant Human Rho GTPase-activating protein 11A protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho GTPase-activating protein 11A (Rho-type GTPase-activating protein 11A), ARHGAP11A, KIAA0013
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGAP11A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.