





KCNIP4 Antibody

Product Code	CSB-PA750882LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6PIL6
Immunogen	Recombinant Human Kv channel-interacting protein 4 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kv channel-interacting protein 4 (KChIP4) (A-type potassium channel modulatory protein 4) (Calsenilin-like protein) (Potassium channel-interacting protein 4), KCNIP4, CALP KCHIP4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNIP4

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA750882LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA750882LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.